SUNBIOSCIENCE

COMPARISON OF AUTOMATED HIGH-CONTENT IMAGE-BASED

PHENOTYPIC SCREENING MODELS FOR DRUG DISCOVERY:



FROM 2D CULTURES TO 3D SPHEROIDS AND ORGANOIDS

Fabien Kuttler (1), Nathalie Brandenberg (2), Sylke Hoehnel (2), Krisztian Homicsko (3), Matthias Lutolf (4) and Gerardo Turcatti (1)

(1) Biomolecular Screening Facility, School of Life Sciences, EPFL, CH-1015 Lausanne, Switzerland (2) SUN Bioscience SA, EPFL Innovation Park, Building D, CH-1015 Lausanne, Switzerland

(3) Department of Oncology, CHUV, CH-1011 Lausanne, Switzerland (4) Laboratory of Stem Cell Bioengineering, School of Life Sciences, EPFL, CH-1015 Lausanne, Switzerland

Translation of promising preclinical candidate compounds into clinical trials) against different models of colorectal disease. The models range live/dead assays (either label-free or with the standard Hoechst / Calcein We show not only how phenotypic analyses can reveal insights into drugs is still too often failing, highlighting the persistent need for from simple 2D cultures of a standard human colorectal carcinoma cell AM / Ethidium Homodimer-1 combination of fluorescent dyes), allowing mechanisms of drug action depending on the model used, but also how line (HCT116 cells) to 3D cultures of spheroids generated from the same effective and relevant disease models used in drug discovery, especially phenotypic image-based analyses in a toxicity screen context. similar and reproducible the results can be between models. However, more subtle differences in effects and efficacy of the tested drugs can for solid tumors. The latest advances in cell-based assay technologies cell line on agarose layers and up to colorectal cancer spheroids and For each model, the obtained images are properly analyzed, segmented patient-derived gastrointestinal organoids using an automated culture be detected only with the state-of-the-art organoid culture derived from and quantified, and hundreds of features are extracted from each of for drug discovery include 3D cultures of patient-derived stem cellbased organoids which are now becoming more standardizable and via stem-cell aggregation in hydrogel-based microcavity arrays (*). them. Standard dimensionality reduction algorithms and multivariate patients, in line with the genetic mutation status of specific oncogenes, automatable for routine use at industrial scale. making this scalable and automatable organoid-culture technology The readouts of screening with the different models are all making use projections are used to integrate the results into simple linear unavoidable in the future of drug discovery through high-content Here we report the direct comparison of screening results of the same classifications of efficiencies of the screened compounds and compared of high-throughput image-based label-free digital holographic set of carefully selected compounds (either FDA-approved or in clinical microscopy or automated fluorescence microscopy and consist of between the different models, as well as confronted to the literature. screening of tissue-level function.



U

Stra

0

()



Unique small molecule library screened

	Scre	en s	et of	80 compounds			s for	for 96-well			plates format		
	1	2	3	4	5	6	7	8	9	10	11	12	
Α	DMSO	Capeci tabine	Caboza ntinib	Ralti trexed	Tofaci tinib	VRT75 2271	Erlo tinib	Olapa rib	Suni tinib	Campto thecin	Tunica mycin G	Gambo gic Acid	
в	DMSO	Cobime tinib	Entrec tinib	Gane tespib	BYL- 719	LEE 011	Bosu tinib	Dasa tinib	Rapa mycin	BAY 43- 9006	Podoph y Ilotoxin	Gambo gic Acid	
С	DMSO	LY2157 299	Tivo zanib	RG 7388	Pacri tinib	Oran tinib	Azacy tidine	Lapa tinib	Afa tinib	Cedi ranib	Taxol	Gambo gic Acid	
D	DMSO	CAL- 101	Lenva tinib	LDE 225	SAR24 5409	Varli tinib	Vande tanib	Dovi tinib	Enza staurin	lmiqui mod	NT 2 Toxin	Gambo gic Acid	
E	DMSO	Dacomi tinib	BLZ945	Birina pant	INCB 024360	Vismo degib	PLX- 4032	Dina ciclib	Fore tinib	Regora fenib	RA- XIII	Gambo gic Acid	
F	DMSO	Bati mastat	Pona tinib	XL- 147	Afure sertib	SAHA	Lenali domide	Tivan tinib	Pexidar tinib	Oxali platin	Colce mid	Gambo gic Acid	
G	DMSO	Apa	Defac	Nera	LY283	Nilo	Gefi	Pozio	Evero	Masa	Vinde	Gambo	





A set of 80 drugs (FDA-approved or in clinical trials) has been selected based on updated clinical data on colorectal cancer treatment. The set is subdivided according to the known target or mechanism of action (MOA) of the molecules into 5 main categories (antimetabolites, toxins / apoptosis inducers, kinase inhibitors, cytokines or immune cells targeting, and Tyrosine Kinase Receptor inhibitors (EFGR, FGF, VEGFR...)

The screen was performed on each model at 1 μ M and 10 μ M final concentration after 2 or 3 days of incubation with the drugs, followed by automated QPI and HCS reading

t-SNE of all screened compounds per model

- t-SNE algorithm was applied on all screening data for all models tested and color-coded according to the MOA of compounds.
- The size of the dots reflects the score of each compounds for the specific model tested, as specified for each graph.





2D culture **HCT116**

HCT116 cells cultured in 2D, and screened as exemplified with OPD DHM). show the reading (through inducing toxic drugs apoptosis i.e. vindesine, colcemid, podophyllotoxin, or specific TK inhibitors, i.e. volasertik or the HSP90-inhibitor ganetespib

Pilaralisib

Cabozantinib 🧧

5-Fluorouracil

Batimastat Pexidartinib abozantinib



Enzastaurin Crizotinib

Spheroids (Agarose) Spheroids Gri3D 2D cells - Ratio Fluo 2D cells – Cell count 2D cells – Confl. (DHM) 2D cells – OPD (DHM)

control

Positiv



Clusterized Heat Map of screening results for the different models tests

Hierarchical clustering / Dendrogram for all tested compounds show overall high similarity in responses of the different tested models

Only few compounds show particular effects on specific models. The main observed difference concerns the subset of compounds targeting the receptor Tyrosine kinase (i.e. EGFR) on the arrayed organoids from CRC patients cultured on Gri3D plates

Receptor Tyr Kinase targeting compounds subset (EGFR mainly, VEGFR...)



Scatterplot matrix of response to drugs for each model

Antimetabolites Kinase inhibitors



Conclusions

However, more subtle differences in The same carefully selected set of effects and efficacy of the tested drugs compounds adainst different models can be detected only with the state-ofcolorectal cancer. ranging from 2D cultures of HCT116 cells the-art organoid culture derived from to **3D** spheroids and patients, in line with cultures the genetic of gastrointestinal organoids of colorectal mutation status of specific oncogenes, cancer patient using an automated 3D making this scalable and automatable culture (Gri3D[®], SUN Bioscience*). organoid-culture technology We show highly similar and reproducible unavoidable in the future of drug results between the models tested. discovery.

show the most intense response with the same set of compounds as in 2D cultures, plus some additional targets i.e. PI3k (alpelisib, voxtalisib)





CRC patient cells cultured in 3D in multiple subwells of a Gri3D plate, forming several dozens homogeneous organoids per well, show the most intense response with the same set of compounds as in 2D cultures, plus some specific targets i.e. EGFR inhibitors (osimertinib...) and other TK inhibitors

belinostat, volasertib, enzastaurin, pacritinib, auranofin, crizotinib, tivantinib or staurosporine EGFR 🔴 ТК 😑 antimetabolite 😑 apoptosis immune 5-Azacytidine Birinapantinib Grantinib RA-XIII CeditaBiertib Vemurafenib Epacadosta Wildlimod Nilotinib Legalidonide Lapatinib Nilotini Turgalidonide Lapatinib Sonidesfectinib Idasanutlin Olapartismodegib Varlitinib Idelalisib Apatinib Pilaralisib Batimastat 🔍 Sunitini Pexidartinib ribociclib Cabozantinib 📍 Everolim Repart vcin



863-874 (2020)