


Unveiling Drug Responses in Liver Spheroids: Multiplexing 3D Cell-Based Assays and Imaging in a Microwell Platform

Ballabio, M.¹, Sirugue, P.², Clapés Cabrer, M.¹, Meyer, M.², Brandenburg, N.², Hoehnel-Ka, S.¹

¹SUN bioscience SA, Lausanne, Switzerland - ²Doppl SA, Lausanne, Switzerland

INTRODUCTION



Yet most 3D cultures are still not ready for screening. Using conventional plates such as ultra-low attachment (ULA) U-bottom poses **challenges** in:



1. Reproducibility
2. Scalability
3. Imaging
4. Automation
5. Handling
6. Loss of microtissues
7. Assay development

To tackle these hurdles, we introduce Gr₃D[®]: an innovative hydrogel microwell 96WP system designed for uniform cell seeding, efficient aggregation, and the generation of individual microtissues in suspension-like conditions. Importantly, the resulting microtissues are strategically positioned within the same focal plane, enabling simultaneous high-resolution imaging. In our study, we cultivated hepatocyte spheroids simultaneously in

Gr₃D[®] and ULA 96WP to compare ease of use, performance, sensitivity, and scalability across multiple assays. We concentrated our efforts on hepatotoxicity. Using Spark Cyto imaging multimode plate reader, we concurrently assessed the viability and enzymatic activities of the 3D microtissues, combining cell imaging and luminescence readouts. This innovative 3D multiplexing approach enabled us to glean crucial insights into the effects of various drugs on CYP activity and microtissue viability, all within a single experiment.

METHODS

Upcyte[®] hepatocyte spheroids are generated in parallel in Gr₃D[®] 96WP 600 μm microwells (SUN bioscience) and in standard ULA 96WP plates (Corning). Spheroids form from a single-cell suspension and are cultured for up to 5 days. The microtissues are exposed to inhibitors of liver enzyme activity for 24 hours, and both CYP (cytochrome P450) activity and microtissue viability are assessed using commercially available kits (Promega). Both plates are evaluated on a Spark[®] Cyto (Tecan, For Research Use Only). All the experiments were carried out in Doppl's Laboratory (Lausanne, Switzerland).



RESULTS




Figure 1) Comparison between Upcyte[®] hepatocyte spheroids grown in parallel A. in Gr₃D[®] 96WP 600 μm microwells or B. in ULA 96WP plates. Both sets of spheroids originate from the same single-cell suspension and are cultured for up to 5 days with a seeding density of 7'500 cells per well. Images captured with a Tecan Spark[®] Cyto - 4x magnification. Scale bars: 500 μm.

Upcyte[®] hepatocyte cells aggregate and form compact spheroids over a 5-day culture period. We seed 7'500 cells per well, resulting in a single spheroid in each ULA well, and 55 microtissues of ~136 cells in each Gr₃D[®] well (Fig. 1). In Gr₃D[®], the total surface area exposed to media and reagents is more than ten-fold greater due to the presence of multiple smaller microtissues (Fig. 2). Consequently, we observe a higher ATP concentration measured from an equivalent total cell number (Fig. 3). Upon exposure to CYP enzyme inhibitors Ketoconazole (CYP3A4) and α-Naphthoflavone (CYP1A2), the viability of microtissues remains unchanged, evidenced by the constant LDH release (Fig. 4). Notably, spheroids in Gr₃D[®] show much higher CYP3A4 basal activity compared to those in ULA. When treated with the well-known inhibitor, Ketoconazole, this activity can be halved in Gr₃D[®]. However, such a reduction remains undetectable in ULA. As for CYP1A2, its basal enzymatic levels are virtually non-detectable; thus, the inhibitory effect of α-Naphthoflavone is not discernible (Fig. 5).




Figure 2) Quantification of spheroid perimeters. A. Images depict the designated region of interest, as analysed with ImageJ. B. Table summarising perimeter values in both Gr₃D[®] and ULA. Scale bars: 500 μm.




Figure 3) Viability of Upcyte[®] hepatocyte spheroids grown either on Gr₃D[®] or ULA 96WP plates (CellTiter-Glo[®]3D). Two-tailed T-test, **P<0.01, ns: non-significant.




Figure 4) Cytotoxicity of Upcyte[®] hepatocyte spheroids grown either on Gr₃D[®] or ULA 96WP plates when treated with CYP inhibitors. A. Brightfield images of spheroids, both treated and untreated, reveal no discernible difference. B. LDH release over time of Upcyte[®] hepatocyte spheroids grown on either Gr₃D[®] or ULA 96WP plates (LDH-Glo[™]).




Figure 5) CYP activity measurements of Upcyte[®] hepatocyte spheroids grown either on Gr₃D[®] or ULA 96WP plates (P450-Glo[™] CYP3A4 and CYP1A2). Left: CYP3A4 activity inhibition upon treatment with Ketoconazole. Right: CYP1A2 activity inhibition after treatment with α-Naphthoflavone. Two-way ANOVA with Fisher's multiple comparisons, *P<0.05, P***<0.001, ns: non-significant.

CONCLUSIONS

- Utilizing a 96-well plate format, Gr₃D[®] 600 μm microwells facilitates the creation of over 50 microtissues per well, allowing for **scalability, uniformity, and robustness**.
- Both **spheroid viability and enzymatic activity detection are improved in Gr₃D[®]**, which positions the platform as an optimal choice for advanced 3D-based drug evaluations.
- The synergistic use of Gr₃D[®] technology with a compatible multimode plate reader, together with multiplexed cell-based assays, enriches our understanding of **liver spheroid responses to drug exposure**.
- Our pioneering approach offers a promising solution to address long-standing challenges in **large-scale 3D cultures and compound evaluation** on physiologically relevant models.

REFERENCES

1. Proctor, WR, et al., 2017, 'Utility of spherical human liver microtissues for prediction of clinical drug-induced liver injury', *Arch Toxicol.*, 91(8):2849-2863, doi: 10.1007/s00204-017-2002-1
2. Brandenburg, N., Hoehnel, S., Kuttler, F., et al., 2020, 'High-throughput automated organoid culture via stem-cell aggregation in microcavity arrays', *Nat Biomed Eng* 4, 863-874, https://doi.org/10.1038/s41551-020-0565-2
3. Tolosa, L., et al., 2016, 'Human Upcyte Hepatocytes: Characterization of the Hepatic Phenotype and Evaluation for Acute and Long-Term Hepatotoxicity Routine Testing', *Toxicological Sciences*, Volume 152, Issue 1, Pages 214-229, https://doi.org/10.1093/toxsci/kfw078

ACKNOWLEDGEMENTS

We thank L. Bruzzone (Tecan) and K. Haupt, K. Oelstrom, V. Pierroz, and A. Fignolé (Promega) for their support. Tecan kindly granted access to the Spark[®] Cyto and Promega provided the assay kits.