





## STANDARD

Gri3D® 96-wellplate standard catalog items are available in 4 microwell sizes (400, 500, 600, 800 µm), and 2 plate bottom options, plastic and imaging bottom.


### Gri3D® STANDARD - PLATE FORMATS

96 wells in 96 wellplate




**Black-walled imaging Gri3D®**  
500 µm microwells (73 per well)  
IBIDI polymer bottom (180 µm)

### Gri3D® STANDARD MICROWELL SIZES AND SAMPLE PLATE



#### Gri3D® sample plate

Available only for standard sizes.  
Contains 1 column filled with hydrogel for each of the standard microwell sizes (400 - 800 µm).



Maximum order size = 3 sample plates.

### Gri3D® - PLATE BOTTOM

Gri3D® 96 wellplate are available in **plastic** (1 mm thick) as well as **imaging bottom** (0.18 mm, IBIDI polymer).

## CUSTOM


### Gri3D® CUSTOM - MICROWELL SIZES

Microwell diameter (in µm)	100	200	300	400	500	600	800	1000	1600	3000
96 wellplate microwell number (per well)	1677	511	211	121	73	55	31	19	7	-
24 wellplate microwell number (per well)	5681	1573	733	379	241	163	102	61	19	7


### Gri3D® CUSTOM - PLATE FORMATS

Gri3D® custom plates can be ordered on 96 wellplate or 24 wellplate formats in a column-basis. One single microwell size can be chosen per plate. Production upon demand with 5% extra fee, minimal order size 5 units per purchase order.


Multiple columns in 96 wellplate




96 wells in 96 wellplate



60 wells (no edges) in 96 wellplate



24 wells in 24 wellplate



### Gri3D® CUSTOM - PLATE BOTTOM

Gri3D® 96 wellplate are available in **plastic** (1 mm thick) as well as **imaging bottom** (0.18 mm, IBIDI polymer).  
Gri3D® 24 wellplate are available in **plastic** (1 mm thick).