


STANDARD


Gri3D® - MICROWELL SIZES




Gri3D® - PLATE FORMATS

Gri3D® standard catalog items are available in 96 wellplate format in 3 configurations.

24 wells in 96 wellplate




96 wells in 96 wellplate



Gri3D® sample plate

Available only for standard sizes. Contains 1 column filled with hydrogel for each of the standard microwell sizes (400 - 800 µm).



Maximum order size = 3 sample plates.

Gri3D® - PLATE BOTTOM

Gri3D® 96 wellplate are available in **plastic** (1 mm thick) as well as **imaging bottom** (0.18 mm, IBIDI polymer).

CUSTOM


Gri3D® CUSTOM - MICROWELL SIZES

Microwell diameter (in µm)	100	200	300	400	500	600	800	1000	1500	3000
96 wellplate microwell number (per well)	1677	511	211	121	73	55	31	19	7	-
24 wellplate microwell number (per well)	5681	1573	733	379	241	163	102	61	19	7


Gri3D® CUSTOM - PLATE FORMATS

Gri3D® custom plates can be ordered on 96 wellplate or 24 wellplate formats in a column-basis. One single microwell size can be chosen per plate. Production upon demand with 5% extra fee.


8 wells in 96 wellplate



96 wells in 96 wellplate




4 wells in 24 wellplate



From 1 column ...
to a full plate

24 wells in 24 wellplate



Maximum 3 24 wellplate per order.

Gri3D® CUSTOM - PLATE BOTTOM

Gri3D® 96 wellplate are available in **plastic** (1 mm thick) as well as **imaging bottom** (0.18 mm, IBIDI polymer).

Gri3D® 24 wellplate are available in **plastic** (1 mm thick).